

中性子散乱実験装置についてもう一度考える —積み木型の中性子散乱実験装置—

北大院工 古坂道弘、佐藤他加志、藤田文行、本間彰、
奥沢康裕、石川直樹、杉田宰、武田晋、宮田環、中西誠一
京大原子炉実験所 杉山正明、日野正裕

In preparation...



中性子散乱実験装置についてもう一度考える —積み木型の中性子散乱実験装置—

2

北大院工 古坂道弘、佐藤他加志、藤田文行、本間彰、
奥沢康裕、石川直樹、杉田宰、武田晋、宮田環、中西誠一
京大原子炉実験所 杉山正明、日野正裕



Conventional SANS instrument

3

Neutron SAS instruments

4

- Typical SANS instruments
- about 10 mmØ sample size



SANS-U@JRR-3




New D11 @ ILL
<http://www.ill.eu/instruments-support/instruments-groups/instruments/d11/news-from-d11/>

Small Pin-hole Time-of-flight SANS

5

- Usually, sample diameter $\approx 10 \text{ mm } \varnothing$ at SANS-U
- 16 m+16m at maximum
- $Q_{\min} \approx 5 \times 10^{-3} \text{ \AA}^{-1}$




Electron linac based pulsed neutron source at Hokkaido University

45MeV Electron Linac @Hokkaido University

- The first generation **compact pulsed neutron source**
 - First beam ≈1975
 - still running...
- 35 MeV, 30 μ A, 50 pps


Hokkaido University Electron Linac



Time averaged intensity

8

Factor 2-3 uncertainties



Tailor made experiment.
example 1.


Small Pin-hole Time-of-flight SANS

Small Pin-hole Time-of-flight SANS

10

- What do you think of this SANS instrument at Hokkaido Univ.
 - Sample size $\approx 2 \text{ mm } \phi$
 - flight path only 1.5m
 - no vacuum tube
- Poor intensity $\approx 1/25???$



$$I(Q) \propto \phi \cdot d\Omega_i \cdot \frac{d\Sigma}{d\Omega} \cdot V_{\text{sample}} \cdot \eta \cdot d\Omega_f$$



Very clean direct beam!

11


- Hokkaido Univ. small pin-hole SANS, 2 mm \varnothing
 - Cd plate drilled with Boric acid



SANS near the direct beam

Preliminary

12




- Cf. Conventional SANS machine:
 - Direct beam stopper covers low angle part! Cf. $\approx 0.004 \text{ Å}^{-1}$ with 8 m
 - SANS B.G. due to collimation

Data analysis

13

- 2mm \varnothing sample SANS is possible even at Hokkaido University.
- TOF necessary
- Direct beam profile is precisely defined
 - Triangular shape direct beam
 - With fine resolution area detector
- How to de-convolute the result?
 - Fit with resolution function convoluted with $S(Q)$?





Ex. 2. Diffractometer

Very low resolution

Diffraction in a steel plate

Preliminary

15




Ex. 3. Intermediate-angle scattering

In between small-angle and wide angle.

Intermediate-angle scattering

17



- $0.05 \text{ \AA}^{-1} < Q < 4 \text{ \AA}^{-1}$
 - with $\lambda=1, 12 \text{ \AA}$
- Ex. Precipitations in steel

What is the tailor-made experiment (TME) beamline

Compact pulsed neutron source is a
best place for TME beamline.

Moderator development


19

- Best place for developing high-performance moderator
 - Coupled moderator concept
 - Solid methane moderator
- Hands-on modification!
 - Low flux means low activation
- J-PARC moderators are developed here.
 - SNS composite moderator tested.



Just imagine some components

20



Hokkaido University
electron linac




Pulsed cold neutron source
solid methane@ \approx 20K

Neutron beam exit

Just imagine some components

21





Tables "optical bench"


Supermirror components

22

ellipsoidal mirror




Bent supermirror



Si perfect crystal bender

23

- High resolution area detector (RPMT)
- Resistive wire type Photomultiplier



Shielding blocks!


24

Tape-engineering!



Detector components

25



- Linear Position Sensitive Detectors

(LPSDs)

8 LPSDs/module \times 6 units



- He-gas:
difficult to purchase?



- High resolution area
detector
(RPMT)
Resistive wire type
Photomultiplier

Resistive wire type PMT +ZnS scintillator

26

- Li (n, α); ZnS(Ag) scintillation
- 3inch, 5inch PMT
 - R2486-04
- Good resolution
 - <1mm
- convenient system for optical device test.


Hirota, Satoh et al.
(RIKEN, KEK, NOP)



mfSANS at Hokkaido Univ.

27

Data acquisition system Detector Sample Ellipsoidal mirror Beam port



Ex. 4. High intensity monochromator

for mfSANS@JRR-3

Commonly used monochromator


29

- Pyrolytic Graphite
 - Mosaic crystal
 - Reflectivity $\approx 80\%$
 - $\Delta\lambda/\lambda \leq 1\%$
 - depending on beam divergence

Not suitable to
SANS, Reflectometer

Large beam divergence

- λ spread
- mosaic




$$\begin{array}{l} \lambda > \lambda_0 \\ \lambda = \lambda_0 \\ \lambda < \lambda_0 \end{array}$$

Collimation effects

- λ spread smaller
- Lower intensity

$$\begin{array}{l} \lambda > \lambda_0 \\ \lambda = \lambda_0 \\ \lambda < \lambda_0 \end{array}$$



PG

Neutron SAS instruments

30

- 1 instrument at the end of beamline.



SANS-U@JRR-3




New D11 @ ILL
<http://www.ill.eu/instruments-support/instruments-groups/instruments/d11/news-from-d11/>

SANS monochromator

31


- Velocity selector is the choice
 - High intensity quasi-monochromator
 - $\Delta\lambda/\lambda \approx 5 \sim 20 \%$
- 1 selector/Guide-end
- Supermirror/monochromator would be other possibility
 - NOP beamline

Extremely
good
performance



Test experiment at H.U. pulsed neutron source


32



Wide λ band realized

33

- λ dependence by time-of-flight measurement
 - Not easy at a reactor.




High intensity monochromator with a beam deflector

New concept!

34

- Improved monochromator being tested





Idea of compact SANS instruments



35

- Many SANS modules
 - along a beamline
- To have long enough machine time anytime,
 - Ex. For protein solution research
- Educate young scientists
- Check your new ideas
- Your own SANS instrument



Tailor made experiment beamline

36

- Test new ideas
 - Very easy to tailor an instrument.
- Best way of learning neutron scattering and instrumentation.
 - To educate young people
 - Teach yourself
- Low cost
- Relaxed resolution to gain intensity
- Tailor made experiment:
 - Ex. just to look at...
 - Intermediate $Q \approx 0.1\text{--}2 \text{\AA}^{-1}$ region
 - a magnetic diffraction peak

Just to review focusing
SANS instrument

Neutron SAS instruments

38

- Need a cold neutron source.
- Only 1 instrument at a guide-end.



SANS-U@JRR-3



New D11 @ ILL
<http://www.ill.eu/instruments-support/instruments-groups/instruments/d11/news-from-d11/>

X-ray SAS instruments are Small...

X-ray SAS instrument

40

- Laboratory sources are available anywhere
- Proper training methods are established
 - Variety of books available
 - You can ask experienced people around.
- If laboratory source is not appropriate, you can always use synchrotron radiation



SAXS

[http://matrix-dev.ansto.gov.au/research/
bragg_institute/facilities/instruments/saxs/
mimi_saxs](http://matrix-dev.ansto.gov.au/research/bragg_institute/facilities/instruments/saxs/mimi_saxs)



ESRF

[http://www.lightsources.org/cms/?
pid=1000103](http://www.lightsources.org/cms/?pid=1000103)



Why not have compact neutron SAS instrument

Idea of compact SANS instruments



42

- To have long enough machine time anytime,
 - Ex. For protein solution research
- Educate young scientists
- Check your new ideas
- Many SANS modules
- Your own SANS instrument



to match with large neutron facilities

43

- Big facilities like
 - ILL, Oak Ridge, Munich, NIST, JRR-3, KAERI, ANSTO, ...
 - ISIS, Luhan center, PSI, SNS, J-PARC



J-PARC
<http://j-parc.jp/ja/layout-j.html>

How to make compact instruments: Focusing is the key.

44

Focusing makes instrument **compact!**
compact \approx low cost

Cf. Focusing gives us **no-gain in intensity!**

Focusing SANS

45

- Toroidal mirror focusing to extend low-Q limit.
 - Moved from Jülich to München
 - $Q=4\times10^{-4}\text{A}^{-1}$




Fig. 1. Toroidal mirror with the image in the detector plane.
B. Alefeld et al./ Physica B 234-236 (1997) 1052-1054

- MgO_2 lens, sextupole lens are available.

Ellipsoidal mirror focusing SANS

Kamada et al. (Hokkaido Univ.)

46

- Ellipsoidal mirror
- 1~10 mm ϕ aperture




図 4.12 Aperture : 1mm ϕ


図 4.17 Aperture : 10mm ϕ

Focusing SANS instrument is Compact!


47

- Focusing

\approx compact



- Conventional point collimation



- Same resolution/intensity

- Angular resolution


$$\approx D/L \approx d/\ell$$

- Intensity:


$$I \propto \phi \cdot d\Omega_i \cdot \frac{d\Sigma}{d\Omega} \cdot V_{sample} \cdot \eta \cdot d\Omega_f$$

Focusing in k space

48



Reciprocal space



Mini-Focusing Small-Angle Neutron Scattering Instrument (mfSANS)

@Hokkaido University

49

Time-of-flight focusing SANS

mfSANS at Hokkaido Univ.

50


Data
acquisition
system

Detector

Sample


Ellipsoidal
mirror

Beam
port



Bovine thighbone, cross section SANS preliminary analysis

51



mfSANS@JRR-3/JAEA



52

Monochromatic neutron
focusing SANS

Prototype focusing SANS@JRR-3


53

- Ellipsoidal mirror
 - 2.5 Q_c supermirror
 - 2.5 m between focal points
 - short radius 20 mm



mfSANS@JRR-3

54




Ni powder 20nm

Preliminary data

55

- $Q_{\min} = 5 \times 10^{-3} \text{ \AA}^{-1}$ using 2mmØ aperture.





wider-angle scattering

Preliminary data




56

- 48 Linear position sensitive detectors at higher angle
 - 1/2 inch dia, 600 mm in length
 - GE made



Micro-phase separated block copolymer


57



Accelerator-driven Laboratory-size neutron-source

58

- RFQ accelerator + DTL \approx 3-11MeV
 - Li or Be target
- Combined with:
 - many mfSANS modules,
 - mini-reflectometers,
 - mini-powder machines??



Summary

59

- Tailor-made neutron scattering instruments beamline
 - To check new ideas
 - Good for self-education
- Small-pinhole SANS demonstrated
- Various new devices under development
 - High performance monochromator
- Very compact SANS instruments under development
 - at Hokkaido university (pulsed neutron source)
 - at JRR-3 reactor of JAEA
 - High-performance SANS modules of your own at large facilities
- With focusing technique
we can build very compact SANS instrument

